PRODUCT PERFORMANCE AT ROOM TEMPERATURE CONDITIONS

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Most of the Eurocell controls and calibrators require 2 to 8 Degrees Celsius storage, typically the temperature range of a refrigerator. This was chosen to extend product shelf life using a controlled environment easily found in most clinical labs worldwide. There was no formal study to identify the ideal temperature range of these products; a temperature study was executed to assess the performance of a product left at room temperature for an extended period of time.

Protocol

Two protocols were used. Protocol 1 used six different product types of whole blood based controls and calibrators and were exposed to 15°C (+/- 2 degrees) for a period of 7 days. Protocol 2 exposed the same six product types of whole blood based controls and calibrators to 20°C (+/- 2 degrees) for a period of 7 days. The pass/fail criteria for both protocols was that all results of testing on the assayed parameters must stay within the assay ranges during the open vial testing performed at the end of the stated shelf life for that product.

Test Results

All values remained with the assay ranges for product stored at 15°C (+/- 2 degrees) for a period of 7 days.

Excluding random errors, all values remained with the assay ranges for product stored at 20°C (+/- 2 degrees) for a period of 7 days with the exception of the mean cell volume (MCV) on the low level of controls. When stored at this temperature, the MCV increases and increased outside of assay range for the low level.

Conclusion

Product that is left outside of the 2-8°C storage temperature for extended periods of time may show an increase in the MCV, depending on the exposure temperature and length of time.

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R & D Systems, Inc.