

R&D Systems CBC-3D™ Assay Table For Sysmex XP-Series • Level Low

R&D systems
a biotechnie brand

EuroCell
Diagnostics

Lot **B0224-L**



05.05.2024



| Param | Unit | Mean | Limit | Barcode | Param | Unit | Mean | Limit | Barcode |
|-------|----------------------|------|-------|---------|--------|----------------------|------|-------|---------|
| WBC | 10e ³ /μl | 1,90 | 0,40 | | LYM | 10e ³ /μl | 1,1 | 0,8 | |
| RBC | 10e ⁶ /μl | 2,35 | 0,16 | | MXD | 10e ³ /μl | 0,1 | 0,1 | |
| HGB | g/dL | 5,9 | 0,4 | | NEUT | 10e ³ /μl | 0,8 | 0,8 | |
| HCT | % | 16,8 | 2,5 | | RDW-CV | % | 12 | 5 | |
| MCV | fl | 71,5 | 6 | | RDW-SD | fl | 35 | 11 | |
| MCH | pg | 25,1 | 2,5 | | MPV | fl | 8,4 | 3,5 | |
| MCHC | g/dl | 35,1 | 4 | | PDW | fl | 10,5 | 10 | |
| PLT | 10e ³ /μl | 60 | 20 | | PCT | % | 0,05 | 0,04 | |
| LYM | % | 57,5 | 42,5 | | P-LCR | % | 17 | 15 | |
| MXD | % | 2 | 2 | | | | | | |
| NEUT | % | 40,5 | 40,5 | | | | | | |

EC REP



Bio-technie® - 19 Rue Louis Delourmel,
CS 49228 Noyal Châtillon Sur Seiche, 35092 Rennes Cedex 9, France

R & D Systems, Inc. - 614 Mc Kinley Place N.E. - Minneapolis,
MN USA 55413

R&D Systems CBC-3D™ Assay Table For Sysmex XP-Series • Level Normal

R&D systems
a biotechnie brand

EuroCell
Diagnostics

Lot **B0224-N**



05.05.2024



| Param | Unit | Mean | Limit | Barcode | Param | Unit | Mean | Limit | Barcode |
|-------|----------------------|------|-------|---------|--------|----------------------|------|-------|---------|
| WBC | 10e ³ /μl | 7,20 | 1,00 | | LYM | 10e ³ /μl | 2,2 | 2,2 | |
| RBC | 10e ⁶ /μl | 4,50 | 0,20 | | MXD | 10e ³ /μl | 0,2 | 0,2 | |
| HGB | g/dL | 13,1 | 0,5 | | NEUT | 10e ³ /μl | 4,8 | 3,4 | |
| HCT | % | 36,2 | 2,7 | | RDW-CV | % | 10 | 5 | |
| MCV | fl | 80,4 | 6 | | RDW-SD | fl | 34 | 12 | |
| MCH | pg | 29,1 | 3 | | MPV | fl | 8,8 | 3,5 | |
| MCHC | g/dl | 36,2 | 4 | | PDW | fl | 10,5 | 10 | |
| PLT | 10e ³ /μl | 270 | 40 | | PCT | % | 0,24 | 0,15 | |
| LYM | % | 31 | 31 | | P-LCR | % | 20,5 | 15 | |
| MXD | % | 2,5 | 2,5 | | | | | | |
| NEUT | % | 66,5 | 33,5 | | | | | | |

EC REP



Bio-technie® - 19 Rue Louis Delourmel,
CS 49228 Noyal Châtillon Sur Seiche, 35092 Rennes Cedex 9, France

R & D Systems, Inc. - 614 Mc Kinley Place N.E. - Minneapolis,
MN USA 55413

R&D Systems CBC-3D™ Assay Table For Sysmex XP-Series • Level High

R&D systems
a biotechnie brand

EuroCell
Diagnostics

Lot **B0224-H**



05.05.2024



| Param | Unit | Mean | Limit | Barcode | Param | Unit | Mean | Limit | Barcode |
|-------|----------------------|-------|-------|---------|--------|----------------------|------|-------|---------|
| WBC | 10e ³ /μl | 18,60 | 2,00 | | LYM | 10e ³ /μl | 2,6 | 2,6 | |
| RBC | 10e ⁶ /μl | 5,70 | 0,25 | | MXD | 10e ³ /μl | 0,6 | 0,6 | |
| HGB | g/dL | 17,8 | 0,7 | | NEUT | 10e ³ /μl | 15,4 | 5,2 | |
| HCT | % | 49,2 | 3 | | RDW-CV | % | 10 | 5 | |
| MCV | fl | 86,3 | 6 | | RDW-SD | fl | 36,5 | 13 | |
| MCH | pg | 31,2 | 3,2 | | MPV | fl | 8,7 | 3,5 | |
| MCHC | g/dl | 36,2 | 4 | | PDW | fl | 10,5 | 10 | |
| PLT | 10e ³ /μl | 615 | 60 | | PCT | % | 0,54 | 0,25 | |
| LYM | % | 14 | 14 | | P-LCR | % | 19,5 | 15 | |
| MXD | % | 3 | 3 | | | | | | |
| NEUT | % | 83 | 17 | | | | | | |

EC REP



Bio-technie® - 19 Rue Louis Delourmel,
CS 49228 Noyal Châtillon Sur Seiche, 35092 Rennes Cedex 9, France

R & D Systems, Inc. - 614 Mc Kinley Place N.E. - Minneapolis,
MN USA 55413